



Lovango S.L. C/ Llull, 63 08005 Barcelona  
tel 933 005 284 fax 934 855 992  
info@lovango.com www.ovan.es

# Jar Test



## ***OVAN multiblade shaker for flocculation.***

- Digital flocculator with 4 or 6 positions controlled by microprocessor.
- Adjustable speed from 20 to 250 rpm with digital indicator, controlled by servo + encoder.
- Smooth and progressive start and stop.
- Digital display shows selected speed and time.
- Speed could be selected precisely every 5 rpm.
- 2 speed values are stored in memory.
- Selection by membrane push-button.
- Digital timer allows a time selection of 0 - 240 min.
- Stainless steel stirring rods of 6 mm Ø with height adjustment.

- Backlight screen that could be placed in 2 positions: lower or rear, without dismantling.
- Suitable for long vessels of 1 L.
- Handles for easy transport.
- Quiet operation. suitable for continuous running.

Dimensions : 80 x 20 x 41 cm. Weight : 15 Kg.

*Ref: 10000-00017 Jar Test*

## **Accessories**

**Injector of coagulants for reagent (6 and 4 places)**

*Ref: 20000-00079*